


Sucrose gradient fractionation

AF Aleksra Filipovska

Updated date: Mar 25, 2021

 An abbreviated version of this protocol was published in Science Advances in Dec 2019

Fidelity of translation initiation is required for coordinated respiratory complex assembly

DOI: 10.1126/sciadv.aay2118

Detailed protocol

1 mg of heart mitochondria, were lysed in 1% n-dodecyl β -D-maltoside (Sigma). Lysates were loaded on 10–30% sucrose gradients and separated by centrifugation overnight at 71, 000 g at 4°C. Gradient fractions were collected as 750 μ l aliquots. RNA was extracted from one third of each fraction by using the TRIzol LS Reagent (Invitrogen), according to the manufacturer's recommendations. The samples were subsequently treated with DNase I and used for cDNA synthesis (Qiagen). The transcript abundance in each fraction was assessed by qRT-PCR. The remainder of each fraction (500 μ l) was precipitated with 72% trichloroacetic acid, resolved by SDS- PAGE and OXPHOS or ribosome-containing fractions were detected by immunoblotting using antibodies specific for individual proteins.

For details on how to pour a sucrose gradient please follow the link below, it is all the same except we use 30% instead of 50% sucrose:

<https://drummondlab.org/protocols/protocol/sucrose-gradient>

How to cite: (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Filipovska, A. (2021). Sucrose gradient fractionation. Bio-protocol Preprint. bio-protocol.org/prep973.
2. Rudler, D. L., Hughes, L. A., Perks, K. L., Richman, T. R., Kuznetsova, I., Ermer, J. A., Abudulai, L. N., Shearwood, A. J., Viola, H. M., Hool, L. C., Siira, S. J., Rackham, O., Filipovska, A. and Filipovska, A. (2019). Fidelity of translation initiation is required for coordinated respiratory complex assembly . Science Advances 5(12). DOI: [10.1126/sciadv.aay2118](https://doi.org/10.1126/sciadv.aay2118)

Copyright: Content may be subjected to copyright.